



PzFlex[®]



Modelling and Fabrication of High Frequency Ultrasound Transducer Arrays for Medical Applications

Robert T. Ssekitoleko

DTC in Medical Devices, University of Strathclyde




Gerry Harvey

Weildlinger Associates, Ltd (PzFlex Europe)

Christine E.M. Démoré, Zhen Qiu, Muhammad. R. Sadiq, Sandy Cochran

Institute of Medical Science and Technology, University of Dundee

UIA 2011: Glasgow, UK, 23-25 May, 2011



Outline

Background

- Ultrasound Systems and Applications
- High Frequency Imaging
- Ultrasonic Device Development

Miniature Array Design

- Design Specification

Fabrication Techniques

- Bonding and Interconnects

Single Crystal Piezoelectric Materials

- Piezoelectric Material Characterisation

Virtual Prototyping

- Preliminary Results
- 15 MHz virtual prototyped transducer


Basic Fabricated Prototype

Summary & Conclusions

Further Work

Ultrasound Systems and Applications


Output to clinicians



Input from technologists

Ultrasound Systems and Applications

Output to clinicians



Input from technologists

High Frequency Imaging

Increasing operational frequency
improves image resolution

- 15 MHz and above give image resolution of less than 0.5 mm

Main applications include:

Intravascular imaging


Small animal imaging

Ophthalmology

Dermatology

BUT


Increasing frequency reduces depth of penetration due to frequency dependent attenuation




Skin layers imaged with 32 MHz transducer. The epidermis (E), Dermis (D) and subcutaneous (S) layers are clear

Transducer on Interventional Tool

Conventional
Ultrasound Probe



Miniature array on a
biopsy needle



Transducer on Interventional Tool

Transducer size reduces with increasing frequency

- 15 MHz array has dimensions small enough to fit in a tool with a diameter of 2 mm

Imaging transducers on interventional tools overcomes the attenuation problem


- Can be placed at the tissue region of interest

Accurate and quick characterisation of tissues with high resolution images


- Miniature arrays enable real-time imaging with electronic scanning of ultrasound beam

Potential for ultrasound *in vivo* pathology


- Accurately positioned tissue biopsies for staging disease




Ultrasonic Device Development



Miniature Array Design



Design Specification




Application:	Imaging of small, deep-seated tumours
Interventional tool:	2 mm diameter biopsy needle
Array Type:	Linear array for electronic scanning Side viewing for rectilinear image
Image resolution:	100 μm or better
Imaging frequency:	15 MHz or higher


Design Specification

Initial Array Specification

- Piezoelectric layer:** 1-3 piezocomposite with single crystal PMN-PT and hard-setting epoxy
- Matching layer:** Alumina loaded epoxy
- Backing layer:** Tungsten loaded epoxy
- Array elements:** 64 elements
100 μm pitch (λ at 15 MHz)
- Array dimensions:** 0.8 mm wide
6.4 mm long
- Cabling:** Flex circuit cabling to fit in core of 2 mm biopsy needle




Basic Array Implementation



- Process diagram is fully defined
- Technical challenges have been outlined and solutions investigated
- Packing technique for the needle orientation is established



Fabrication Techniques




Array Fabrication

1-3 Piezocomposite is made from fragile **PMN-PT piezocrystals**


Composites made through **dice and fill** to enhance properties

Lapped and polished to a known thickness for the desired frequency

Matching and backing layers are cast on to **avoid bondlines**




1-3 piezocomposite




Interconnect

Imaging array



Bonding and Interconnects



Fragile PMN-PT materials need **low temperature and pressure** for bonding electrical interconnects


Conductive silver-loaded epoxy cures at room temperature

Electromagnetic anisotropic UV curable epoxy has a good potential


Photolithographic patterning of flex-circuit with **dry film photoresist** and copper etching

190 μm track width and 110 μm separation achieved as first prototype

Technique shows promising results for prototypes with finer pitch tracks



Single Crystal Piezoelectric Materials



PMN-PT Single crystal

Benefits

Improved properties compared to conventional ceramics

- For imaging
 - Improved coupling coefficient
 - High permittivity

Drawbacks

Fragile in fabrication

Expensive


Difficult to characterise

- Properties vary even if from the same batch


Makes modelling more challenging

Parameter	Symbol	PZT-5H	PMN-PT
Relative permittivity at constant strain	ϵ^{s33}	1470	3026
Electromechanical coupling coefficient	k_t	0.51	0.62
Curie temperature	T_c	195 °C	150 °C

Piezoelectric Material Characterisation




45° rotated k31 plates k31 plates




IEEE standard characterisation technique

- Multiple specimen geometries of plates and bars to isolate multiple resonant modes
- Electrical impedance spectroscopy of specimens analysed with PRAP (TASI Technical Software Inc., Kingston, Canada)
- Full set of elastic, piezoelectric and dielectric properties can be extracted




Uniformity of Material Properties



PMN-PT material from the same supplier

Varied material properties




(a)

For example relative permittivity at constant strain varies by **~80%**

Z.Qui et al., "Characterisation of piezocrystals for practical configurations with temperature and pressure-dependent electrical impedance spectroscopy." IEEE Trans. Ultrason., Ferroelect. and Freq. Control, in press 2011.

Transducer Design with Virtual Prototyping




Virtual Prototyping

Finite element analysis (FEA) allows **realistic and cost-effective** analysis of **device performance**


Virtual prototyping enables **rapid testing** of ranges of design parameters

PZFlex, a time domain FEA software utilised for the design

However, FEA requires **accurate material properties** for a realistic mode




2D Composite with matching




3D Composite

Virtual Prototyping

- **First model:**
 - Monolithic block of piezoelectric material
 - Material properties are modified to match simulation to experimentally obtained electrical impedance spectra
- **Second model:**
 - 1-3 piezocomposite block
 - Check electrical impedance from simulation and experiment match for confidence in full model
- **Third model:**
 - Full array with 1-3 composite, matching and backing layers
 - Can also model bond lines or electrode thicknesses to determine effects of fabrication processes
 - Used to optimise design



2D Block




2D Composite

Preliminary Results: Model 1

The model was validated with a 2D plate resonating at low frequency

PTZ-5H (CTS 3203HD)


10x10x1 mm




PMN-PT Model was more challenging to match to the experiment matches with experimental test

2D plate model was used to extract properties


10x10x1 mm plate



Preliminary Results: Model 1



Supplier A



Supplier B

- It is more challenging to get a one model fits all for piezocrystal



Preliminary Results: Model 2

1-3 PMN-PT/Epoxy Composite


3d unit cell model

5 MHz resonant frequency

Some properties were increased by up to 80% to obtain a reasonable match




Fitted curves on linear scale




Fitted curve on log scale

15 MHz virtual prototyped transducer



- Active layer thickness is 80 μm
- Matching layer thickness is 50 μm
- Backing thickness is 800 μm



- 36% piezocrystal volume fraction in the composite has higher impedance magnitude at the electrical resonant frequency than the 49%.

Basic Fabricated Prototype

Single element transducer integrated into small tube


15 MHz transducer with matching and backing layers

Microfabrication techniques such as lapping and micro-dicing were used

Transducer element: 1 mm x 5 mm

Metal tube diameter: 2 mm

Coaxial cable diameter: 0.34 mm





Summary & Conclusions

High frequency transducers and arrays for **high resolution imaging** have significant clinical potential

- Miniature arrays enable integration with interventional tools
- Design of device and probe can be optimised for specific applications

Fabrication of a prototype **array integrated into a biopsy needle** is in progress


Comprehensive **characterisation of single crystal piezoelectric materials** is required for more accurate results in virtual prototyping

Virtual prototyping is used to save time and expense in optimising the design of high frequency transducers



Further Work

- Prototype testing
 - Impedance analysis
 - Pulse-echo test
 - Imaging
- More material characterisation
- More device fabrication working through the process diagram
 - Fabrication of small pitch flexi-circuit
- Ex-vivo tissue testing



Acknowledgements



Prof. George Corner and Dr. Elaine Henry
(Ninewells Hospital, Dundee, UK)

MicroSystems Engineering Centre
(Heriot-Watt University, Edinburgh, UK)



Dr. Jeff Bamber and Nigel Bush
(Joint Department of Physics,
Institute for Cancer Research, London)

Loadpoint Ltd. (Swindon, UK)



AFM Ltd. (Birmingham, UK)

Logitech Ltd. (Glasgow, UK)

